



EC DECLARATION OF CONFORMITY

(For safety of personal protective equipment 89/686/EEC)

Reference : DSD151201 / 01-12-2015

This is certify that the following designated products

- Item description : **FINGER TAB**
- Model : 622600 , 622800 , 622400 , 622810 , 622430
622500 , 622700 , 622300

Conforms to the following Directives and Standards:

89/686/EEC (Safety)

- **EN 388 : 1994**

TEST REPORT NO: TAK-013788

This Declaration of Conformity is responsibility of the manufacturer

Manufacturer Name : DOOSUNG CORPORATION

**Address : 29 42 Sorae-ro, Daegot-Myeon, Gimpo-si,
gyeonggi-do, Korea.**

I, the undersigned, hereby declare that the equipment specified above,
conforms to the above Directive(s) and Standard(s):

DOO SUNG CORPORATION

D. Y. Lee. President



Signature _____

Date : December. 01, 2015

Name : Dall Young, Lee



YOUR PARTNER FOR THE BEST QUALITY

TEST REPORT

98, Goyukwon-ro, Gwacheon-si, Gyeonggi-do, 427-716, Rep. of KOREA TEL 82-2-2164-0011 FAX 82-2-2634-1008

Report No : TAK-013788

Receipt Date : Aug.27.2015

Client : DALL YOUNG, LEE

Test Completion Date : Sep.04.2015

DOOSUNG CORPORATION

928 Song Mar-Ri, Daegot-myeon, Kimpo-city, Kyungki-Do, Korea.

Sample : FINGER TAB (622600/622800/622400/622810/622430/622500/622700/622300)

TEST RESULTS

TEST ITEM	UNIT	SAMPLE	RESULT	TEST METHOD
Cr(VI)	mg/kg	-	Not Detected	EPA 3060A

- Method Detection Limit -

Cr(VI) : 0.5 mg/kg

* USAGE : SUBMITTED PURPOSE FOR

NOTE : 1. The test results on this test report are only limited to the samples and sample names provided by the customer and KTR does not guarantee the quality of all products of the customer, and you can confirm the authenticity of the test report online (www.ktr.or.kr) or by using the QR code.
2. This test report shall not be used for public relation, advertisement, lawsuit and any other purposes outside the scope of its defined usage.
3. Only original copy (including certified copy) of the test report is valid - electronic copy (hard and/or soft) is for your reference only.

Hong Chang-ki

Prepared by Hong Chang-ki
E-mail: environchem@ktr.or.kr

Kim Sun-il

Reviewed by Kim Sun-il
Technical Manager
Tel: 1577-0091(ARS ①-④)

Sep.04.2015

Korea Testing & Research Institute

President

Choi Hyeon-ki



QR Code to verify genuineness

1 of Total 1 Page(s)

622300**MIDAS FINGER TAB**

MIDAS FINGER TAB is a full-featured finger tab with genuine cordovan leather face. It features adjustable finger spacer and shelf for precise tuning. The shiny anodized base plate is designed to fit the palm of the hand improving finger stability at full draw.

SPECIFICATION

► MATERIAL	Aluminum plate / Cordovan leather face
► SIZES	SS / S / M / L
► TYPE	RH / LH

**622600****SMART FINGER TAB**

SMART FINGER TAB is high quality leather finger tab. It features anodized aluminum plate and replaceable cordovan leather face. Includes shelf and finger spacer.

SPECIFICATION

► MATERIAL	Aluminum plate / Cordovan leather face
► SIZES	S / M / L
► TYPE	RH / LH

**622400****CARTEL FINGER TAB**

CARTEL FINGER TAB is made of high quality cordovan leather face with anodized aluminum plate backing to provide a maximum comfort and good finger protection while shooting.

SPECIFICATION

► MATERIAL	Aluminum plate / Cordovan leather face
► SIZES	M / L
► TYPE	RH / LH



622500**CR-306 FINGER TAB**

High-performance tab design improves release feel. CR-306 features genuine cordovan leather face with metal plate backing for perfect finger protection and stable shooting.

SPECIFICATION

▶ MATERIAL	Aluminum plate / Cordovan leather face
▶ SIZES	S / M / L
▶ TYPE	RH / LH

**622700****PRO FINGER TAB****SPECIFICATION**

▶ MATERIAL	Synthetic leather
▶ SIZES	M / L
▶ TYPE	RH / LH

**622430****PRO-I FINGER TAB****SPECIFICATION**

▶ MATERIAL	Aluminum plate / Synthetic leather
▶ SIZES	M / L
▶ TYPE	RH / LH

**622810****PRO-III FINGER TAB****SPECIFICATION**

▶ MATERIAL	Synthetic leather
▶ SIZES	S / M / L
▶ TYPE	RH / LH

**622800****PRO-II FINGER TAB****SPECIFICATION**

▶ MATERIAL	Synthetic leather
▶ SIZES	S / M / L
▶ TYPE	RH / LH

